

PALFIQUE OMNICHROMA

Every shade, Every class, One choice


Technique Guide

Anterior layering technique
By: Mustafa Fahmi
& Omnicroma and Blocker technique





OMNICHROMA BLOCKER

Anterior




Posterior



Leading the Way
in Dental Innovation and
Unmatched Experience

PALFIQUE LX5



OPACITY



OPA2 OA2 A2
 0.75 0.65 0.55



ESTELITE COLOR



Enamel

Dentine

Palatal shield
 Palfique LX5 WE
 Omnichroma
 Palfique LX5 CE

Dentine
 Palfique LX5 OPA2

Body
 Palfique LX5 OA1
 Palfique LX5 OA2
 Palfique LX5 OA3
**according to the shade of the tooth*

Enamel
 Palfique LX5 A1
 Palfique LX5 A2
 Palfique LX5 A3
 Palfique LX5 A3.5
 Palfique LX5 B1
 Palfique LX5 B2
 Palfique LX5 C1
 Palfique LX5 C2
 Palfique LX5 CE
 Palfique LX5 WE
 Omnichroma
**according to the shade of the tooth*

Colors
 Estelite color Clear
 Estelite color White
 Estelite color Blue
 Estelite color Pink Opaque
 Estelite color Gray
 Estelite color LCO
 Estelite color MCO
**according to the shade of the tooth*

Dentine

Enamel

